

sempermed®
supreme plus

DISTRIBUTION PAR



STERIGENE
INGÉNIERIE DES PROCÉDÉS PROPRES & STÉRILES

SUPREME PLUS MORE GRIP TO YOUR WORLD

EXCELLENT GRIP – perfect instrument guidance
ABSOLUTELY SKIN-FRIENDLY – the Semperit Post Production Process
OPTIMUM SAFETY – ideal tactile properties



MORE GRIP TO YOUR WORLD

The Sempermed® Supreme Plus provides you with more grip in every situation while combining maximum safety, ideal tactical properties and perfect instrument guidance.

■ ABSOLUTELY SKIN-FRIENDLY

The specially developed, automatic finishing process contains intensive and multiple washing steps. Skin-irritating substances resulting from the production process and the natural latex protein concentration are thus reduced to the technically lowest level.

■ THE LINING SYSTEM

The special net-like lining structure optimises the contact surface between the film and the skin, thus considerably reducing the friction resistance when putting on the glove. The result: The Sempermed® Supreme Plus facilitates easy donning of the gloves.

■ DOUBLE DONNING

For your safety and easy donning use Sempermed® Supreme as under glove and Sempermed® Supreme Plus as outer glove. Recommendation: If you wear different sizes, put on the larger pair first. Or use gloves of the same size.



TECHNICAL DATA

sempermed® supreme plus	
TYPE	Sterile surgical glove for single use, powder-free, with synthetic lining
MATERIAL	Natural rubber latex (NR)
COLOUR	Natural white
GLOVE SHAPE	Fully anatomical with rolled rim
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ min. 270 mm 7, 7½ and 8 min. 280 mm 8½ and 9 min. 285 mm
WALL THICKNESS measured in single layer in the palm area	0.19–0.23 mm
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead ≥ 9 N)	≥ 12 N
DURABILITY	3 years in original packaging if stored as per DIN 7716, ISO 2230
STERILISATION	Radiation STERILE R
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Directive 89/686/EEC	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical fibre-free peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs
ARTICLE NUMBERS	SPTSUPLUS. Tailles: 5.5 / 6.0 / 6.5 / ... / 9.0



CONTACT US!



2 rue André Citroën
95130 FRANCONVILLE
Tél. : +33 (0)1 34 44 23 23
Fax : +33 (0)1 34 44 23 24
www.sterigene.com

¹⁾Please note the Sempermed Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock.**